



BIOPTRON Compact III Light Therapy device is designed for use at home, in hospitals and in therapeutic centers. The device is easy to handle for anybody, and in order to provide comfort it is equipped with a timer that gives a signal every two minutes. The device is packed in a handy plastic case. The filter diameter of the BIOPTRON Compact III is approx. 4 cm.

TECHNICAL DATA

Filter diameter approx.	4 cm
Power supply	100-230 V~, 50/60 Hz
Power consumption	56 VA
Rated power of halogen	20 W
Protective class	Class II, IP 20
Device type	Type BF 
Weight without Tripod	0.5 kg
Ambient temperatures	
- Operation	+10°C to +40°C
- Storage	-30°C to +45°C
Wavelength	480 - 3400 nm
Degree of polarization	>95% (590 - 1550 nm)
Specific power density	av. 40 mW/cm ²
Light energy per minute	av. 2.4 J/cm ²
CE labeling	

In the United States, contact Easy Figure International for information at www.EasyFigure.com

In Canada, contact Easy Figure International for information at www.EasyFigure-Zepter.com