

TECHNICAL DATA-Bioptron 2

Filter diameter approx.	15 cm
Power supply	100-240 V~, 50/60 Hz
Power consumption	1.4 - 1.0 A
Fuse	T2A / 250 V
Rated power of halogen	90 W
Protective class	Class I, IP 20
Weight without stand	4.3 kg
Ambient temperatures	
- Operation	+10°C to +40°C
- Storage	-30°C to +45°C
Wavelength	480 - 3400 nm
Degree of polarization	>95% (590 - 1550 nm)
Specific power density	av. 40 mW/cm ²
Light energy per minute	av. 2.4 J/cm ²
CE labeling	

