

TECHNICAL DATA – BIOPTRON PRO 1

Wavelength	480 - 3400 nm
Degree of polarization	>95% (590 - 1550 nm)
Specific power density	<i>av. 40 mW/cm²</i>
Light energy per minute	<i>av. 2.4 J/cm²</i>
Light Intensity	min. 10.000 lux
CE labeling	
Weight without stand	3.4 kg
Weight with stand	7.8 kg
Energy consumption in sleep mode	0.5 W
Power Supply	100-240 V~ 50/60 Hz
Power consumption	90 VA

