


TECHNICAL DATA, YOUTHRON/MEDALL

| | |
|--------------------------------|---|
| Filter diameter approx. | 5 cm |
| Power supply | 100-240 V~, 50/60 Hz |
| Power consumption | 0,29-0,12 A |
| Rated power of halogen | 20 W |
| Protective class | Class II <input type="checkbox"/> , IP 20 |
| Weight without stand | 0.5 kg |
| Ambient temperatures | |
| - Operation | +10°C to +30°C |
| - Storage | 0°C to +40°C |
| Wavelength | 480 - 3400 nm |
| Degree of polarization | >95% (590 - 1550 nm) |
| Specific power density | av. 40 mW/cm ² |
| Light energy per minute | av. 2.4 J/cm ² |
| CE labelling |  |

