

- *Effective*
- *Safe*
- *Non-invasive*
- *Painless*
- *Swiss Made*



BIOPTRON YOUTHTRON

*PROVEN PATENTED LIGHT TECHNOLOGY
for health and beauty, for all skin types*

BIOPTRON YouTHron

PAG-970



The BIOPTRON Color Light Set provides well-being for health-conscious people. Using the seven colors of the light spectrum, BIOPTRON Color Light aims to balance and enhance your body's energy centres and to stimulate your body's own healing process. BIOPTRON Color Light Therapy is a result of a wide experience in the therapeutic field of color and light. It works mainly in two ways:

1. BIOINFORMATION on PRESSURE POINTS - color applied to your skin in a focused form, using bioinformation technology to activate your cells via pressure points.
2. STIMULATION OF CHAKRA ZONES, on the different body and mind centres. The BIOPTRON Color Light Set consists of 7 colored glass filters and 3 bioinformed therapeutic cosmetic products which work together to enhance any color therapy session. The color filters are hand-made in France from high-quality glass, ensuring optimum refraction of light and brilliant color intensity. Besides being handmade with the utmost skill and care, the glass and cosmetic products are enriched with bioinformation.



Technical data

Light features	<ul style="list-style-type: none"> • Polarized • Polychromatic UV-free • Incoherent • Low energy
Power supply	100–240 V~ 50/60Hz
Power consumption	0,29-0,12 A
Protective class	II □ IP 20
Weight	0,5 kg 2,1kg
Wavelength	480–3400 nm
< 0.3 W in sleep mode	
Degree of polarization	>95% 590–1550nm
Specific power density	max. 40mW/cm ²
Light energy per minute	max. 2,4J/cm ²
Light intensity	~10000 lux
CE labelling	
According to	IEC/EN60335-1 EN/IEC 61000-3-2/-3

